
LDNB2317
Compressor
Nebulizer System

INSTRUCTIONS FOR USE


L'Docile Enterprise GmbH Germany

Mouthpiece


Small Mask

Atomizer


Big Mask

Air Tube


Air Openings

Tube Connector

Power Switch

Filter & Cover

LDNB2317


Compressor Nebulizer System

● COMPONENTS

1. Power switch
2. Air openings
3. Air tube
4. Atomizer
5. Mouth piece
6. Adult mask
7. Child mask
8. Tube connector
9. Filter and cover

● WORK PRINCIPLE

This unit is compact, simple to operate and convenient to carry. It uses a non-oil lubricating, single-cylinder piston pump to create air pressure. The air flow created by the air compressor comes out of the jet outlet and then enters into the nebulizer kit. The main air flow into and out of the air compressor unit is illustrated in the figure below.

Class II Device (double safety insulation).

Type B Device (device with specific protection against electrical hazards).

● INTRODUCTION

The device is an AC powered air Nebulizer Compressor System intended to provide a source of compressed air for medical purposes for use in home healthcare. It is to be used with the original kit nebulizer to produce aerosol particles of medication for respiratory therapy for both children and adults.

Important: Read these instructions carefully before using this device.

● INSTRUCTIONS FOR USE

1. After unpacking the device, check the compressor and accessories for visible damage or defects, such as cracks in the plastic housing which could expose electrical components. Contact customer service to report any damaged product.
2. Before using the device, proceed with the cleaning operations as described in the Cleaning and Disinfection of the accessories section.
3. Open the atomizer by turning the top counter clockwise.
4. Put the prescribed quantity of medication into the medicine cup. Maximum fill volume is 8 ml.
5. Close the atomizer by holding the bottom and turning the top part clockwise, making sure that the unit is thoroughly closed.
6. Connect one end of the air tube to the atomizer, and the other end to the air outlet on the device.
7. Connect the angled mouthpiece to the atomizer.
8. If using the mask, remove the angled mouthpiece and connect the mask directly to the atomizer.
9. Plug the compressor into the electrical outlet, making sure the electrical rating shown on the rating plate on the bottom of the unit corresponds to your main voltage.
10. To start the treatment, switch the device ON by pressing the ON/OFF switch onto the I position.
11. When the treatment has been completed, switch the unit OFF by setting the switch position to O position and disconnect the plug from the electrical outlet.
12. Wash the atomizer and accessories as described in the Cleaning and Disinfection of the accessories section.

⚠ CAUTION

1. Use the device only for its intended purpose as described in these instructions, for use as a compressor nebulizer system for aerosol therapy. Follow the instructions of your doctor. Any other form of use constitutes improper use and may be dangerous. The manufacturer cannot be held liable for any damage caused by improper, incorrect or unreasonable use, or if the unit is connected to electrical installations not complying with safety regulations.
2. Keep this manual for future reference.
3. Do not operate this unit in the presence of any anesthetic, oxygen or inflammable mixtures.
4. In case of failure and/or malfunction, read the Troubleshooting section. Do not open the compressor housing as this is not a serviceable product.
5. The proper function of the device can be compromised if original replacement parts e.g., filters are not used.
6. Never submerge the unit in water as the device is not protected against water penetration.
7. Never touch the unit when your hands are wet or moist.
8. Do not leave the unit outdoors.
9. Place on a stable and horizontal surface when operating.
10. Make sure the air vent openings are not obstructed.
11. Do not allow unsupervised children or physically challenged individuals to use the unit.
12. Do not turn the unit OFF by simply pulling out the plug from the electrical outlet.
13. Do not use any adapter, simple or multiple, and/or extension cable.
14. Do not leave the unit plugged in when not in use. Unplug the device from the electrical outlet when not in operation.
15. Follow the manufacturer's instructions for installing the device. The manufacturer cannot be held responsible for any damage caused by incorrect installation.
16. Device not suitable for anesthesia or pulmonary ventilation system.
17. The power supply cord cannot be replaced by the user.

18. The power supply cord should always be fully unwound in order to prevent dangerous overheating.
19. Before performing any cleaning operation or filter replacement, disconnect the device from the main supply by unplugging it or by turning the main supply switch OFF.
20. If you decide not to use the device, it is recommended to dispose of it according to the local regulations in force.
21. Remember to use this device only with medications prescribed by your doctor, and use the accessories recommended for the therapy.
22. Verify that the medications prescribed by your doctor can be used with common compressor nebulizer systems for aerosol therapy.
23. Life cycle of the device is 5 years.

● **CLEANING AND MAINTENANCE**

a. Cleaning of the compressor

Use a clean dry cloth and mild detergents to clean the outside unit.

ATTENTION: During cleaning, make sure the internal parts of the unit are not in contact with liquids and the power cord is not connected.

b. Cleaning and Disinfection of the accessories

After you have disconnected the air hose, disconnect the mouthpiece.

Open the cap and empty the medicine container.

Wash the atomizer kit with running water or soak in warm water for 15 minutes.

For more thorough, hygienic cleaning, add some vinegar in the water.

NOTE: Do not boil the atomizer kit. Doing so can damage it.

Dry the atomizer kit completely before storing it.

c. Replacement of the atomizer

The atomizer should be replaced after a long period of inactivity, or if it is obstructed by dry medication, dust, etc. With normal use, your atomizer should last at least 6 months

d. Replacement of the filter

The filter is to be replaced after approximately 30 hours of use, or when it appears dirty. Open the filter holder to replace the filter with a new one. Use original filters only. Do not use the device without the filter.

● TROUBLESHOOTING GUIDE

Symptom	Probable Cause	Solution
Too much noise	No filter in the unit	Put a filter in the unit
The Nebulizer does not run	Main compressor switch is off	Turn the unit on
	Atomizer was not cleaned after its previous use	Clean the unit
	Air hose is bent	Remove any bend or kink in the air hose
	The filter is clogged	Replace the filter
	No medicine is left	Add the appropriate amount of medicine prescribed by your physician to the medicine container
	The cap is not attached properly	Reattach the cap
Water drops are formed in the air hose	The medicine was over-full	Adjust the amount of medicine in the medicine container to the appropriate level, connect the air hose to the air compressor, and turn the unit on
	The nebulizer kit was washed but not dried	Open and close the nebulizer kit repeatedly until the drops disappear

If the device does not work properly after attempting these steps, contact Customer Service.

Maintenance and repairs

Never open the device. There are no serviceable parts within the unit. The compressor requires no lubrication or maintenance.

● DEVICE SYMBOLS

I ON (power)

⚡ Type B

□ Class II

O OFF (power)



~ Alternating current


⌘ Separate collection

⚠ ATTENTION: Read the instructions carefully

SPECIFICATIONS

Air Compressor Nebulizer Machine Features	
Power:	220 or 120V AC, 50/60Hz
Power consumption:	Below 65W
Nebulizer flow:	5 L
Compressor flow:	8-10 L
Medicine container capacity:	2-8 ml
Particle size:	0.3 to 6 microns
MMAD:	1.663um
Average nebulization rate:	0.2ml/min
Remains amount:	0.85ml
Noise level:	below 65dBA
Compressor pressure range:	10 to 43.5 Psi
	70 to 300 Kpa
	0.7 ~3.0 bar
Operation pressure range:	11.5 to 16 Psi
	80 to 110 Kpa
	0.8 ~ 1.1 bar
Operating temperature range:	10°C to 40°C (50°F to 104°F)
Operation humidity:	10 to 95% RH
Storage temperature range:	-20°C to 70°C (-4°F to 158°F)
Storage humidity range:	10 to 95% RH
Dimension (L x W x H):	210 x 140 x 95 mm
Product weight:	1.42 kg
Standard accessories:	Compressor Nebulizer Machine
	Mask (large)
	Mask (small)
	Filters x 5
	Air tube
	Medicine container 8 ml
Mouthpiece	

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Shanghai International Holding Corp. GmbH (Europe)